

Materials and Equipment

- Cell sample, e.g. Blood sample (RBC-depleted), PBMCs or murine splenocytes
- Pro5[®] Recombinant MHC Pentamer conjugated to biotin.
- *Either:*
 - Streptavidin magnetic beads to detect the biotin label of the Pentamer, suitable for use with magnetic columns (e.g. Streptavidin Microbeads, Miltenyi Biotec #130-048-102)
 - Magnetic stand e.g. MACS Multistand, Miltenyi Biotec #130-042-303
 - Separation unit e.g. MiniMACS separation unit, Miltenyi Biotec #130-042-102
 - Magnetic columns e.g. MS Columns, Miltenyi Biotec #130-042-201
- *Or:*
 - Streptavidin magnetic beads to detect the biotin label of the Pentamer, suitable for use with an external magnet (e.g. Dynabeads[®] M-280 Streptavidin, Invitrogen #112-05D; Lodestars[™] 2.7 Streptavidin, Polymer Laboratories #6727-1001)
 - Magnetic tube holder
- 96-well cell culture plates, U-bottom
- Complete IMDM (IMDM medium containing 5% fetal calf serum, 10% AB0 serum (for human cells) or 5% murine serum (for mouse cells) and 50 IU/ml IL-2 + penicillin/streptomycin)
- Wash buffer (0.1% sodium azide, 0.1% BSA in PBS) - remove air prior to use
- Benchtop centrifuge with swing-out rotor and appropriate carriers
- Hemacytometer and microscope for cell counting
- Humidified CO₂ incubator, heated to 37°C

Standard Procedures (may be scaled up or down according to requirements and may need adjustment depending on the particular system used)

Column-based isolation of antigen-specific cells using beads from Miltenyi Biotec

- 1. Stain 5-10 × 10⁶ cells with 5 tests biotin-labeled Pentamer in 200 ml wash buffer for 15 minutes at room temperature.**
- 2. Wash samples with 8 ml wash buffer, spin 400 × g for 5-10 minutes in a chilled centrifuge at 4°C, and discard the supernatant.**
- 3. Resuspend cells in 90 ml wash buffer and add 10 ml streptavidin magnetic beads per tube**(or as directed by the manufacturer).
- 4. Incubate for 20 minutes in the refrigerator.**
- 5. Wash cell-bead complexes with 8 ml wash buffer, spin 400 × g for 5-10 minutes in a chilled centrifuge at 4°C and resuspend in 500 ml wash buffer (de-gassed).**
- 6. Meanwhile, wash a column suitable for positive-selection with 500 ml wash buffer (de-gassed) and place on magnetic stand.**
- 7. Load cell-bead complexes onto the column.** Antigen-specific T cells labeled with Pentamer-bead complexes will be retained on the column (positive fraction).
- 8. Collect the negative fraction that elutes from the column, including 3 washes of 500 ml each.**
- 9. Remove column from magnet and flush out the positive fraction by adding 1 ml wash buffer onto the column and applying a plunger (provided with the column).**

10. To obtain a purer antigen-specific cell population, the positive cell fraction may be passed over a second column.

Tube-based isolation of antigen-specific cells using beads from Invitrogen or Polymer Laboratories

1. Stain $5-10 \times 10^6$ cells with 5 tests biotin-labeled Pentamer in 200 ml wash buffer for 15 minutes at room temperature.
2. Meanwhile, calculate the volume of streptavidin beads required for incubation with the cells. Wash this volume of beads 3 times with 3 ml wash buffer, using a magnetic tube holder, then resuspend in the original volume. 5 beads per cell is recommended in the first instance.
3. Wash cells with 8 ml wash buffer, spin $400 \times g$ for 5-10 minutes in a chilled centrifuge at 4°C and discard the supernatant.
4. Resuspend cells in 200 ml wash buffer.
5. Add the pre-washed streptavidin beads.
6. Incubate for 20 minutes at room temperature, mixing every 5 minutes.
7. Bring the volume in the tube up to 2 ml with wash buffer then place in a magnetic tube holder.
8. Leave for 3-5 minutes. If desired, supernatant can be retained for flow cytometric analysis to confirm removal of antigen-specific cells, otherwise discard supernatant.
9. Wash the fraction containing cell-bead complexes (attracted to wall of tube) 3 times with 500 ml wash buffer and discard supernatant.

Cell culture / expansion

1. Spin the isolated cells at $400 \times g$ for 5-10 minutes in a chilled centrifuge and discard the supernatant.
2. Resuspend cells in 100 ml complete IMDM and transfer to one well of a U-bottom 96-well tissue culture plate. Expansion is best if 1×10^5 irradiated autologous feeder cells (negative fraction from bead isolation) are also added to each well.
3. Place the culture plate at 37°C in a humidified CO_2 incubator for desired incubation period.
4. After 8-10 days of expansion the cells should be counted and/or analyzed by flow cytometry to determine the extent of expansion. For flow cytometric analysis, cells should be re-stained with Pentamer for 10 minutes at room temperature, then washed and incubated with anti-CD8 antibody and Fluorescent-labeled Streptavidin for 30 minutes at 4°C .
5. If further expansion is desired, re-stimulation with fresh irradiated autologous feeder cells should be carried out at day 10 - 12 of culture.
6. With careful culture technique and good re-stimulation every 10 - 12 days, the cells may be cultured for 1 - 2 months.

n.b. Cellular stimulation generally results in down-regulation of the T cell receptor. For this reason, Pentamer staining should be carried out on cells that have been rested (i.e. no re-stimulation) for 7 - 10 days. This will ensure a good level of staining.